

# Subthalamic nucleus deep brain stimulation with a multiple independent constant current-controlled device in Parkinson's disease (INTREPID): a multicentre, double-blind, randomised, sham-controlled study

Jerrold Vitek (1), Roshini Jain (2), Lilly Chen (3), Alexander Tröster (4), Lauren Schrock (1), Paul House (6), Monique Giroux (7), Adam Hebb (8), Sierra Farris (9), Donald Whiting (10), Timothy Leichter (11), Jill Ostrem (12), Marta San Luciano (12), Nicholas Galifianakis (12), Leo Verhagen Metman (13), Sepehr Sani (13), Jessica Karl (14), Mustafa Siddiqui (15), Stephen Tatter (15), Ihtsham ul Haq (16), Andre Machado (17), Michal Gostkowski (17), Michele Tagliati (18), Adam Mamelak (19), Michael Okun (20), Kelly Foote (21), Guillermo Moguel-Cobos (22), Francisco Ponce (23), Rajesh Pahwa (24), Jules Nazzarro (25), Cathrin Buetefisch (26), Robert Gross (27), Corneliu Luca (28), Jonathan Jagid (29), Gonzalo Revuelta (30), Istvan Takacs (31), Michael Pourfar (32), Alon Mogilner (33), Andrew Duker (34), George Mandybur (35), Joshua Rosenow (36), Scott Cooper (37), Michael Park (38), Suketu Khandhar (39), Mark Sedrak (40), Fenna Phibbs (41), Julie Pilitsis (42), Ryan Uitti (43), Philip Starr (5)

1. University of Minnesota School of Medicine , Minneapolis, USA 2. Boston Scientific Corporation, Valencia, USA 3. Boston Scientific Corporation, Valencia , USA 4. Barrow Neurological Institute , Phoenix, USA 5. University of California, San Francisco (UCSF), San Francisco, USA 6. Neurosurgical Associates, LLC, Murray, USA 7. Esai INC., Woodcliff Lake, USA 8. Kaiser Permanente of Colorado, Denver, USA 9. Boston Scientific , Seattle, USA 10. Allegheny General Hospital , Pittsburgh, USA 11. Allegheny General Hospital , Pittsburgh, USA 12. University of California, San Francisco (UCSF), San Francisco, USA 13. Rush University Medical Center , Chicago, USA 14. Rush University Medical Center , Chicago, USA 15. Wake Forest University School of Medicine , Winston-Salem, USA 16. Wake Forest University School of Medicine , Winston-Salem, USA 17. Cleveland Clinic, Cleveland, USA 18. Cedars-Sinai Medical Center , Los Angeles, USA 19. Cedars-Sinai Medical Center , Los Angeles , USA 20. University of Florida , Gainesville, USA 21. University of Florida , Gainesville, USA 22. Barrow Neurological Institute , Phoenix, USA 23. Barrow Neurological Institute , Phoenix, USA 24. University of Kansas Medical Center , Kansas City, USA 25. University of Kansas Medical Center , Kansas City, USA 26. Emory University School of Medicine , Atlanta, USA 27. Emory University School of Medicine , Atlanta, USA 28. University of Miami School of Medicine , Miami, USA 29. University of Miami School of Medicine , Miami, USA 30. Medical University of South Carolina , Charleston, USA 31. Medical University of South Carolina , Charleston, USA 32. New York University Medical Center , New York City , USA 33. New York University Medical Center , New York City, USA 34. University of Cincinnati, Cincinnati, USA 35. University of Cincinnati, Cincinnati, USA 36. Northwestern University School of Medicine , Chicago, USA 37. University of Minnesota School of Medicine , Minneapolis, USA 38. University of Minnesota School of Medicine , Minneapolis, USA 39. Kaiser Permanente of California, Sacramento, USA 40. Kaiser Permanente of California, Redwood City , USA 41. Vanderbilt University Medical Center , Nashville, USA 42. Albany Medical Center, Albany, USA 43. Mayo Clinic, Jacksonville, USA

## PUBLISHED

**Lancet Neurology**

## STUDY TYPE

**Multi-center, prospective,  
double-blind, randomized (3:1)  
controlled trial  
with sham control**

## DEVICE

**Vercise™  
Deep Brain Stimulation System**

This summary is created by Boston Scientific and is intended to consolidate the publication for educational use only. All figures created by Boston Scientific based on the published data. One or more authors report funding from Boston Scientific as a financial disclosure. Boston Scientific was the primary sponsor of this study.

NM-719603-AA

Vitek JL, Jain R, Chen L, et al. Subthalamic nucleus deep brain stimulation with a multiple independent constant current-controlled device in Parkinson's disease (INTREPID): a multicentre, double-blind, randomised, sham-controlled study. *Lancet Neurol.* 2020;19(6):491-501. doi:10.1016/S1474-4422(20)30108-3

BACKGROUND

The INTREPID Study assessed improvement in motor function and quality of life in PD patients following bilateral subthalamic nucleus (STN) DBS using the Vercise™ Deep Brain Stimulation System with multiple independent current sources that allows for selective activation of individual contacts on the DBS lead there by permitting a defined distribution of applied current.

METHODS

Study Design:

- Multi-center, prospective, double-blind, Randomized (3:1) with sham control

Principal Investigators:

- Dr. Philip Starr (University of California, San Francisco)
- Dr. Jerrold Vitek (University of Minnesota)

Primary Endpoint:


- Mean difference between active and control groups in ON time w/o troublesome dyskinesia, with no increase in antiparkinsonian medications (LED), from post-implant baseline to 12 weeks post-randomization

Sites & Subjects:

- N = 160 randomized subjects
- 23 Sites

Sites & Subjects:

- Up to 5 years post-randomization



BASELINE DEMOGRAPHICS AND CLINICAL CHARACTERISTICS

CLINICAL CHARACTERISTICS/ DEMOGRAPHICS AT SCREENING	
Age (years) - Mean (SD) n	59.9 (7.95) 160
Gender – Male (n %)	72.5% (116/160)
PARKINSON’S DISEASE RELATED SYMPTOMS	
UPDRS III Scores (meds OFF) – Mean (SD) n	43.4 (9.60) 153
UPDRS III Scores (meds ON) – Mean (SD) n	18.5 (8.26) 157
Disease Duration (years) - Mean (SD) n	10.1 (3.61) 160
PD diary: OFF time ON time with troublesome dyskinesias ON time without dyskinesias ON time with non-troublesome dyskinesias Asleep	 6.91 ± 2.99 hours 4.35 ± 2.63 hours 4.65 ± 2.67 hours 3.65 ± 1.90 hours 7.20 ± 1.47 hours

KEY INCLUSION/EXCLUSION CRITERIA:

Key Inclusion Criteria

- Diagnosis of bilateral idiopathic PD with ≥ 5 years of motor symptoms.
- Modified H&Y≥ 2; UPDRS-III score of ≥ 30 (meds off); at least 33% improvement in UPDRS III following meds
- Greater than or equal to 6 hours of poor motor function (OFF time plus ON time with troublesome dyskinesias) per day as assessed by PD diary
- An appropriate candidate for surgical procedures required for bilateral STN DBS

Key Exclusion Criteria

- Any intracranial abnormality or medical condition that contraindicates DBS surgery
- Have any significant psychiatric condition likely to compromise the subject’s ability to comply with requirements of the study protocol

RESULTS

PRIMARY ENDPOINT


ON Time Without Troublesome Dyskinesias

At the end of the blinded phase (12-weeks):

- A 3.74 ± 4.79 hour improvement in ON time was reported in the active group compared to a 0.72 ± 3.56 hour improvement in the control group compared to post-implant activation


At the end of the blinded phase (12-weeks):

- The difference in mean change from the baseline visit (post-implant) to 3 months post-randomization in increased ON time between the active and control groups was 3.03 ± 4.52 (P < 0.001)




UPDRS III SCORES (STIM ON/MEDS OFF)

51%  
Improvement in UPDRS III at 1 year (n = 160)




IMPROVEMENT IN ON TIME AT 1 YEAR




IMPROVEMENT IN ON TIME  
6.0 hours

- At 1 year, a 6 ± 3.8 hour increase in ON time without troublesome dyskinesias was reported compared with Screening (P < 0.001)
- Higher improvement than other similar studies in the literature

34%  
Improvement in Quality of Life at 1 year



41%  
of programs used more than  
a single contact at one year



Safety Profile

Safety Profile is comparable to Literature reports

Surgical complications across all patients who underwent implantation (n = 196)  
[whether or not they went on to be randomized in the study]

EVENT	NUMBER OF EVENTS (SUBJECTS)	RATE
Infection occurring in first 6 months of surgery associated with partial/total hardware removal (SAE)	7 (7)	7/196 (3.6%)
Symptomatic Perioperative Intracranial Hemorrhage (SAE) occurring during surgery or postsurgical hospitalization	4 (4)	4/196 (2.0%)
Symptomatic Peri Lead Edema (SAE)	6 (6)	6/196 (3.1%)
Return to OR due to electronic malfunction	4 (4)	4/196 (2.0%)
Return to OR due to lead breakage/lead migration	0 (0)	0/196 (0%)
Return to OR due to hardware erosion occurring > 6 mos. after implantation	2 (2)	2/196 (1.0%)

- Deaths: 7 deaths (unrelated to the study-device and/or implant procedure)
- No unanticipated adverse events
- No lead breakage/lead fractures

Rate = Number of subjects with event/total number of implanted subjects\*100

CONCLUSIONS

- The study successfully met the primary endpoint and several secondary endpoints based on outcomes reported during the 12-week blinded period.<sup>1</sup>
- At 1 year, UPDRS III scores (stim on/meds off) improved 51% and ON time without troublesome dyskinesias (PD-diary) increased 6 ± 3.8 hours compared with screening (p < 0.001), respectively.<sup>1</sup>
- Overall Improvement in quality of life, medication reduction and high satisfaction with therapy was maintained
- The overall safety profile of the DBS System was comparable to other published reports

1. Vitek JL, Jain R, Chen L, et al. Subthalamic nucleus deep brain stimulation with a multiple independent constant current-controlled device in Parkinson's disease (INTREPID): a multicentre, double-blind, randomised, sham-controlled study. Lancet Neurol. 2020;19(6):491-501. doi:10.1016/S1474-4422(20)30108-3

Results from different clinical investigations are not directly comparable. Information provided for educational purposes only.

**Indications for Use:** The Vercise™ Deep Brain Stimulation (DBS) Systems are indicated for use in bilateral stimulation of the subthalamic nucleus (STN) as an adjunctive therapy in reducing some of the symptoms of moderate to advanced levodopa-responsive Parkinson's disease (PD) that are not adequately controlled with medication.

**Contraindications, warnings, precautions, side effects.** The Vercise DBS Systems, or any of its components, are contraindicated for: Diathermy as either a treatment for a medical condition or as part of a surgical procedure, Electroconvulsive Therapy (ECT) and Transcranial Magnetic Stimulation (TMS) as the safety of these therapies in patients implanted with the Vercise DBS Systems have not been established, Magnetic Resonance Imaging (MRI), patients who are unable to operate the system, patients who are poor surgical candidates or who experience unsuccessful test stimulation. Refer to the Instructions for Use provided with the Vercise DBS Systems or [BostonScientific.com](http://BostonScientific.com) for potential adverse effects, warnings, and precautions prior to using this product.

All trademarks are the property of their respective owners.

**Caution:** U.S. Federal law restricts this device to sale by or on the order of a physician.

**Boston  
Scientific**  
Advancing science for life™

25155 Rye Canyon Loop  
Valencia, CA 91355 USA

Copyright ©2020 by  
Boston Scientific Corporation  
or its affiliates. All rights reserved.

NM-719603-AA