


BETTOCCHI Hysteroscopes

with Working Channel for 5 Fr. Operating Instruments

The BETTOCCHI hysteroscopes enable the physician to perform both a purely diagnostic procedure with visual examination of the uterine cavity, but also give him the possibility to perform minor surgical interventions.

The BETTOCCHI hysteroscopes have been specially developed for this process. This instrument's small diameter and anatomically adapted outer shape is optimally designed to meet the requirements of outpatient usage.

- Outer diameter of 4 mm only
- Integrated 5 Fr. working channel
- Huge variety of mechanical or bipolar instruments available
- Recommended for hysteroscopic sterilization using the ESSURE® procedure



- Set 1: BETTOCCHI **Hysteroscope, diameter 4 mm**, based on 2 mm telescope, including:
- 26008 BA HOPKINS® II **Forward Oblique Telescope 30°**, diameter 2 mm, length 26 cm, autoclavable
 - 26152 BI BETTOCCHI **Operating Sheath**, diameter 3.6 mm, with channel for 5 Fr. semi-rigid operating instruments, for use with CF operating sheath 26152 BO
 - 26152 BO BETTOCCHI **Continuous-Flow Operating Sheath**, diameter 4 mm for use with operating sheath 26152 BI
 - 26159 UHW **Biopsy and Grasping Forceps**, semi-rigid, double-action jaws, 5 Fr., length 34 cm
 - 26159 SHW **Scissors**, semi-rigid, pointed, single-action jaws, 5 Fr., length 34 cm

EndoFax
Endo @ Mail







KARL STORZ
Endoscopy-America, Inc.
2151 E. Grand Ave.
El Segundo, CA 90245

Phone: (800) 421-0837
Fax: (800) 321-1304
E-Mail: info@karlstorz.com
Web: www.karlstorz.com

- Set 2:** BETTOCCHI **Hysteroscope diameter 5 mm**, based on 2.9 mm telescope, including:
- 26120 BA HOPKINS® II **Forward Oblique Telescope 30°**, diameter 2.9 mm, length 30 cm, autoclavable
 - 26153 BI BETTOCCHI **Operating Sheath**, diameter 4.3 mm, with channel for 5 Fr. semi-rigid operating instruments, for use with CF operating sheath 26153 BO
 - 26153 BO BETTOCCHI **Continuous-Flow Operating Sheath**, diameter 5 mm, for use with operating sheath 26153 BI
 - 26159 UHW **Biopsy and Grasping Forceps**, semi-rigid, double-action jaws, 5 Fr., length 34 cm
 - 26159 SHW **Scissors**, semi-rigid, pointed, single-action jaws, 5 Fr., length 34 cm

When ordering one of the above-mentioned sets, you will be granted an **additional discount of 25% off**.

Optional Mechanical Operating Instruments for use with BETTOCCHI Hysteroscopes

- | | | |
|---|-----------|---|
|  | 26159 H | HESSELING Tenaculum Grasping Forceps , semi-rigid, double-action jaws, 5 Fr., length 34 cm |
|  | 26159 EHW | Scissors , semi-rigid, blunt, single-action jaws, 5 Fr., length 34 cm |
|  | 26159 DHW | Punch , semi-rigid, through-cutting, single-action jaws, 5 Fr., length 34 cm |
|  | 26159 BHW | Biopsy Spoon Forceps , semi-rigid, double-action jaws, 5 Fr., length 34 cm |
|  | 26159 M | BETTOCCHI Myoma Fixation Instrument , semi-rigid, 5 Fr., length 34 cm |
|  | 26159 G | BETTOCCHI Palpation Probe , semi-rigid, 5 Fr., length 34 cm |

EndoFax
Endo @ Mail

KARL STORZ
Endoscopy-America, Inc.
2151 E. Grand Ave.
El Segundo, CA 90245

Phone: (800) 421-0837

Fax: (800) 321-1304

E-Mail: info@karlstorz.com

Web: www.karlstorz.com

**Optional Bipolar Operating Instruments for use with
BETTOCCHI Hysteroscopes**



26159 BE **Bipolar Vaporization Electrode**,
semi-rigid, 5 Fr., length 36 cm



26159 GC **GORDTS/CAMPO Bipolar Ball
Electrode**, semi-rigid, 5 Fr.,
length 36 cm



26158 BE **Bipolar Vaporisation Electrode**,
5 Fr., 90° angled, with chamfer in handle

When ordering the above-mentioned optional items in addition to a set, you
will be granted an **additional discount of 25% off**.



EndoFax
Endo @ Mail

KARL STORZ
Endoscopy-America, Inc.
2151 E. Grand Ave.
El Segundo, CA 90245

Phone: (800) 421-0837
Fax: (800) 321-1304
E-Mail: info@karlstorz.com
Web: www.karlstorz.com

For further information, please also have a look at our current
Manual of Hysteroscopy - Diagnosis and Surgery
Manual of Permanent Hysteroscopic Intratubal Sterilization
EndoWorld GYN 24-3-E/04-2009