

# LAPAROSCOPY

**CLICK<sup>®</sup>line**  
*Hand Instruments*

## *Freedom of Choice*



**STORZ**  
KARL STORZ — ENDOSKOPE

COAGULATE

CUT

**ROTATING KNOB**  
Controlled 360° rotation

**FLUSHING PORT**  
Allows for optimal cleaning  
of the instrument

**ONE-TOUCH  
RELEASE BUTTON**  
Quick and simple disassembly

**45° CAUTERY PIN**  
Reduces electrocautery cord  
interference

**INTERCHANGEABLE TIPS**  
Cost-effective replacement and versatility

**INSULATED SHAFT**  
Space-age polymer resists cuts, tears  
and arcing

**IMPROVED HANDLE DESIGN**  
One handle fits all instrument sizes while  
reducing hand fatigue

**CLICK<sup>®</sup>line**

ROTATE

- Precision craftsmanship combines a solid feel with superior distal tip response
- Design your own instrumentation from over 7,000 possible combinations
- Modular components for convenience and cost-effectiveness
- Interchangeable design allows standardization for multiple surgical applications
- 1-year manufacturer's defect warranty

GRASP

**STORZ**  
KARL STORZ—ENDOSCOPE

ARTICULATE



# CLICK & go OUTER TUBES

## INSULATED



<b>33200</b>	5mm outer diameter, 30cm working length with Luer-lock connector and plug
<b>33300</b>	5mm outer diameter, 36cm working length with Luer-lock connector and plug
<b>33400</b>	5mm outer diameter, 43cm working length with Luer-lock connector and plug
<b>33500</b>	10mm outer diameter, 36cm working length with Luer-lock connector and plug
<b>33600</b>	10mm outer diameter, 43cm working length (for use with 5mm, 43cm insert)

## NON-INSULATED, METAL



<b>33200 M</b>	5mm outer diameter, 30cm working length with Luer-lock connector and plug
<b>33300 M</b>	5mm outer diameter, 36cm working length with Luer-lock connector and plug
<b>33400 M</b>	5mm outer diameter, 43cm working length with Luer-lock connector and plug
<b>33500 M</b>	10mm outer diameter, 36cm working length with Luer-lock connector and plug

# CLICK & go STERILIZATION



The modular design of Clickline makes it fast and easy to thoroughly clean and sterilize the components. In situations where it is preferred to maintain the instrument in the fully assembled configuration, Clickline instruments have been validated for pre-vacuum steam sterilization with the exception of the articulating sets. For specific specifications contact your KARL STORZ representative or Customer Service at (800)-421-0837.

# CLICK & go HANDLES

**33121** Plastic, insulated



**33125** Metal, insulated



**33122** Plastic, insulated, with Manhes swing-away ratchet



**33126** Plastic, insulated, with disengageable ratchet



**33123** Plastic, insulated, with hemostat ratchet



**33125 P** Metal, insulated, non-rotating



**33131** Metal



**33132** Metal with Manhes swing-away ratchet



**33133** Metal with hemostat ratchet



**33141** Metal with finger trigger and disengageable ratchet



**33137** Metal Y-Handle with ratchet



**30131** Metal mini for precision control





# CLICK & go HANDLES



**33137 F** Metal Y-Handle with spring-action



**30132** Metal mini with disengageable ratchet for precision control



**33138** Metal, straight



**33131 P** Metal, non-rotating



**33140** Metal, angulated with ratchet



**33151** Plastic Handle without ratchet, with comfort grip



**33152** Plastic Handle with Manhes ratchet, with comfort grip



**33153** Plastic Handle with hemostat style ratchet, with comfort grip




**33156** Plastic Handle w/disengageable ratchet, with comfort grip





**33161** Metal Handle without ratchet, with comfort grip



# CLICK & go ACCESSORIES



**33120** Rotation wheel,  
package of 5, reusable



**33131 F** Spring for 33131  
handle to provide  
positive action



**26002 M** High frequency monopolar  
cable w/4mm banana  
plug for monopolar HF  
electrosurgical generators,  
reusable, 300cm



**26004 M** High frequency monopolar  
cable w/4mm banana  
plug for monopolar HF  
electrosurgical generators,  
models Aesculap, Berchtold,  
Martin and Siemens,  
reusable, 300cm



**26005 M** High frequency monopolar  
cable w/5mm banana  
plug for monopolar HF  
electrosurgical generators,  
models Karl Storz and Erbe,  
reusable, 300cm



**26006 M** High frequency monopolar  
cable w/8mm banana  
plug for monopolar HF  
electrosurgical generators,  
models Valleylab and Bovie,  
reusable, 300cm



**27717 D** Sterilization tray,  
36cm



**39312 B** Sterilization  
tray, 43cm



**27717 DHS** Tray insert to  
accommodate  
mini and pediatric  
instruments (need  
to purchase  
2 per tray)















**27650 G/5** Cleaning brush  
for sheath, 50cm,  
package of 5















**27652/3** Cleaning brush for  
jaws, package of 3





























TIPS	TIP IDENTIFIER	AVAILABLE LENGTHS						
		2mm X 20cm	3mm X 20cm	3mm X 30cm	3.5mm X 30cm	5mm X 30cm	5mm X 36cm	5mm X 43cm
	Reddick-Olsen heavy dissector & grasper							
UL					⊙	⊙	⊙	⊙
	Reddick-Olsen heavy dissector & grasper							
ULS		⊙						
	Atraumatic grasping forceps							
AF						⊙	⊙	⊙
	Reddick-Olsen dissecting & grasping forceps							
KJ						⊙	⊙	⊙
	Reddick-Olsen dissecting & grasping forceps							
KJS		⊙						
	Dolphin-nose dissecting & grasping forceps							
DN						⊙	⊙	⊙
	Dolphin-nose dissecting & grasping forceps, small							
MP						⊙	⊙	⊙
	Conus grasper, atraumatic, hollow jaws							
CG						⊙	⊙	⊙
	Atraumatic dissecting & grasping wedge shaped							
SD						⊙	⊙	⊙
	Kelly dissecting & grasping forceps							
MD						⊙	⊙	⊙
	Kelly dissecting & grasping forceps							
MDS		⊙				⊙		
	Maryland dissecting & grasping forceps							
ML						⊙	⊙	⊙

TIPS	TIP IDENTIFIER	AVAILABLE LENGTHS						
		2mm X 20cm	3mm X 20cm	3mm X 30cm	3.5mm X 30cm	5mm X 30cm	5mm X 36cm	5mm X 43cm
	Maryland dissecting & grasping forceps							
MLS		⊙						
	Maryland dissecting & grasping forceps, long							
MLL								⊙
	Atraumatic dissecting & grasping forceps, 40mm jaws							
PR								⊙
	Manhes atraumatic dissector & grasper							
MF						⊙	⊙	⊙
	Straight atraumatic dissecting forceps							
MH						⊙	⊙	⊙
	Straight atraumatic dissecting forceps							
MHS		⊙						
	Manhes straight tip dissector, 15mm							
MU						⊙	⊙	⊙
	Manhes conical blunt tip dissector, 10mm							
MV						⊙	⊙	⊙
	Manhes conical blunt tip dissector, 15mm							
MW						⊙	⊙	⊙
	Manhes duck-bill blunt tip dissector, 15mm							
MR						⊙	⊙	⊙
	Manhes duck-bill blunt grasper, 10mm							
MS						⊙	⊙	⊙
	Manhes duck-bill blunt grasper, 15mm							
MT						⊙	⊙	⊙





































# CLICK & go TIPS













TIPS		TIP IDENTIFIER	AVAILABLE LENGTHS						
			2mm X 20cm	3mm X 20cm	3mm X 30cm	3.5mm X 30cm	5mm X 30cm	5mm X 36cm	5mm X 43cm
	Fine-tooth fenestrated dissector & grasper								
ON									
	Fine-tooth fenestrated dissector & grasper								
ONS									
	Right angle dissecting forcep, standard								
R									
	Right angle dissecting forcep, short								
RS									
	Right angle dissecting forcep, large								
RG									
	Manhes micro forceps with partial grasping, short								
MKM									
	Manhes micro forceps with partial grasping, long								
MLM									
	Manhes micro forceps with full grasping, long								
MNM									
	Alain-Grousseau pyloric spreader								
VDS									
	Tan pyloric spreader and grasper								
VRS									
	Double spoon atraumatic grasper with teeth								
D									
	Alligator jaw grasping forcep								
HM									

TIPS		TIP IDENTIFIER	AVAILABLE LENGTHS						
			2mm X 20cm	3mm X 20cm	3mm X 30cm	3.5mm X 30cm	5mm X 30cm	5mm X 36cm	5mm X 43cm
	Fundus grasper, atraumatic								
FG									
	Fundus grasper								
FGS									
	Wavy atraumatic grasper, hollow jaws								
LF									
	Wavy atraumatic grasper, hollow jaws, double-action								
LFD									
	Wavy grasper with distal teeth								
WT									
	Wavy grasper with distal teeth, double-action								
WTD									
	Schneider lymph node grasper, atraumatic								
SN									
	Babcock grasping forceps								
A									
	Babcock grasping forceps								
AS									
	Babcock fenestrated grasper, short, multi teeth								
BA									
	Babcock fenestrated grasper, long, multi teeth								
BB									
	Babcock grasping forceps								
BC									




















































TIPS	TIP IDENTIFIER	AVAILABLE LENGTHS							
		2mm X 20cm	3mm X 20cm	3mm X 30cm	3.5mm X 30cm	5mm X 30cm	5mm X 36cm	5mm X 43cm	10mm X 36cm
DISSECTORS/GRASPERS									
	Babcock paddle grasper, large square, long	BL							
	Babcock paddle grasper, small round	BLS							
	Duval grasping forceps	DU							
	DeBakey atraumatic grasping forceps	DY							
	DeBakey atraumatic grasping forceps, double-action	DYD							
	Narrow, curved dissecting & grasping forceps, serrated	DK							
	Large organ grasping forceps, 70mm long	SG							
	Bowel clamp with two rows of atraumatic teeth	GF							
	Bowel grasping forceps	CB							
	Bowel grasper, fenestrated	C							
	Grasping forceps	TF							
	Croce-Olmi atraumatic grasping forceps, curved	CC							

TIPS	TIP IDENTIFIER	AVAILABLE LENGTHS							
		2mm X 20cm	3mm X 20cm	3mm X 30cm	3.5mm X 30cm	5mm X 30cm	5mm X 36cm	5mm X 43cm	10mm X 36cm
DISSECTORS/GRASPERS									
	Fenestrated double action grasping forceps	K				●	●	●	●
	Fenestrated double action grasping forceps	KS		●					
	Vancaille fenestrated adhesion forceps	VK					●	●	●
	Atraumatic dissecting & grasping forcep	DF					●	●	●
	Vancaille meso-ovarian forceps	VA					●	●	●
	Mangeshikar serrated grasping forceps	AG					●	●	●
	Mangeshikar grasping forceps, without serration	AR					●	●	●
	Swolin fallopian tube grasping forceps	SW					●	●	●
	Vancaille oviduct forceps	VT					●	●	●
	Nuwayhid grasper	NS		●					
	Fenestrated grasping forceps, short	TH							●
	Matkowitz grasping forceps with distal teeth	KW					●	●	●







































# CLICK & go TIPS







































TIPS	TIP IDENTIFIER	AVAILABLE LENGTHS							
		2mm X 20cm	3mm X 20cm	3mm X 30cm	3.5mm X 30cm	5mm X 30cm	5mm X 36cm	5mm X 43cm	10mm X 36cm
DISSECTORS/GRASPERS									
	Manhes micro grasping forceps	MEM							
	Manhes grasper with multiple teeth	M							
	Manhes cobra-jaw with 1 x 2 teeth	MA							
	Manhes cobra-jaws with 1 x 2 teeth	MAS							
	Grasping forceps, 2 x 4 teeth	G							
	Claw forceps, 2 x 3 teeth, long	FS							
	Claw forceps, 2 x 3 teeth	FM							
	Claw forceps, 2 x 3 teeth, short	FK							
	Manhes tiger-jaws, 2 x 4 teeth	MG							
	Manhes tiger-jaws, 2 x 4 teeth	MGS							
	Sawalhe tissue grasping forceps with bayonet	UN							
	Manhes multiple teeth atraumatic grasper	ME							

TIPS		TIP IDENTIFIER		AVAILABLE LENGTHS							
				2mm X 20cm	3mm X 20cm	3mm X 30cm	3.5mm X 30cm	5mm X 30cm	5mm X 36cm	5mm X 43cm	10mm X 36cm
DISSECTORS/GRASPERS											
	Manhes multiple teeth atraumatic grasper	MES									
	Berci multifunctional grasper, 1 x 1 teeth	J									
	Remorgida grasping forceps, 1 x 2 teeth	RA									
	Remorgida grasping forceps, insulated jaws	RM									
	Sawalhe tenaculum forceps, long	UM									
	Sawalhe tenaculum forceps, short	UO									
	Shawki cutting & grasping forceps	SH									
	Koh ultramicro I grasping forceps, straight	FA									
	Koh ultramicro II grasping forceps, straight	FB									
	Koh ultramicro I grasping forceps, curved left	FL									
	Koh ultramicro I grasping forceps, curved right	FR									
	Spoon forceps	GS									





TIPS		TIP IDENTIFIER	AVAILABLE LENGTHS							
			2mm X 20cm	3mm X 20cm	3mm X 30cm	3.5mm X 30cm	5mm X 30cm	5mm X 36cm	5mm X 43cm	10mm X 36cm
DISSECTORS/GRASPERS										
	Padded graspers for anvils & choledochoscopes	PG								
	Anvil grasper	AV								
SCISSORS										
	Metzenbaum curved scissors, 12mm	MS								
	Metzenbaum curved scissors	MSS								
	Curved scissors with pointed jaws	MSP								
	Curved serrated scissors, 17mm	MA								
	Curved serrated scissors, conical	MW								
	Curved serrated scissors, conical	MWS								
	Straight scissors	MD								
	Manhes serrated scissors, heavy blade	MT								
	Scissors, heavy blades w/cylindrical grinding	WS								

TIPS	TIP IDENTIFIER	AVAILABLE LENGTHS							
		2mm X 20cm	3mm X 20cm	3mm X 30cm	3.5mm X 30cm	5mm X 30cm	5mm X 36cm	5mm X 43cm	10mm X 36cm
SCISSORS									
	Shawki cutting & grasping forceps	SH							
	Hook scissors	EH							
	Micro hook scissors	EHS							
	Micro hook scissors	ES							
	Hook scissors, jaws not crossing	EK							
	Self sharpening micro scissors, long	SF							
	Self sharpening micro scissors, short	SK							
	Micro dissecting scissors	S							
	Koh ultramicro suture scissors, sharp	SB							
	Koh ultramicro suture scissors, sharp, short	SBS							
	Koh ultramicro suture scissors, blunt	SA							
	Micro dissecting scissors, curved	SL							



# CLICK & go TIPS

		AVAILABLE LENGTHS						
		TIP IDENTIFIER	2mm X 20cm	3mm X 20cm	3mm X 30cm	3.5mm X 30cm	5mm X 30cm	5mm X 36cm
TIPS								
<b>BIOPSY &amp; PUNCH</b>								
	Blakesley dissecting and biopsy forceps							
	DB							
	DBS							
	Biopsy punch forceps, through-cutting							
	DH							
	DHS							

		AVAILABLE LENGTHS						
		TIP IDENTIFIER	2mm X 20cm	3mm X 20cm	3mm X 30cm	3.5mm X 30cm	5mm X 30cm	5mm X 36cm
TIPS								
<b>BIOPSY &amp; PUNCH</b>								
	Biopsy punch forceps with two teeth							
	DZ							
	Manhes biopsy forceps							
	MB							
	Frangenheim biopsy punch, through-cutting							
	DS							

# CLICK & go BARIATRICS

Meet the growing demand for surgical treatment of morbid obesity with a comprehensive line of bariatric instrumentation from KARL STORZ. Advances in minimally invasive laparoscopy now yield similar promising patient outcomes as open surgery - but with shorter hospital stays, better cosmesis, and faster recovery.

Whether you are considering performing gastric bypass, LAP-BAND, or general surgery on an obese patient, KARL STORZ has the latest, high-precision instrumentation to broaden and extend your clinical capabilities.

- 43cm extended length instruments
- 42cm extended length Hopkins® II laparoscopes
- Anvil graspers
- Atraumatic graspers
- EndoTIP® trocars ideal for obese patients







BAYONET



CURVED DOWN



CURVED UP

# CLICK & go CURVED

TIPS	TIP I.D.	STYLE	5mm X 30cm	5mm X 36cm
<b>GRASPERS</b>				
	Cuschiери grasper with multiple teeth, long			
	<b>BBO</b>	Curved up	<input type="radio"/>	<input type="radio"/>
	<b>BBU</b>	Curved down	<input type="radio"/>	<input type="radio"/>
	Cuschiери grasper with multiple teeth, short			
	<b>BAB</b>	Bayonet	<input type="radio"/>	<input type="radio"/>
	<b>BAO</b>	Curved up	<input type="radio"/>	<input type="radio"/>
	<b>BAU</b>	Curved down	<input type="radio"/>	<input type="radio"/>
	<b>MB</b>	Bayonet	<input type="radio"/>	<input type="radio"/>
	<b>MO</b>	Curved up	<input type="radio"/>	<input type="radio"/>
	<b>MU</b>	Curved down	<input type="radio"/>	<input type="radio"/>
	Cuschiери grasper with multiple teeth, 4.8mm wide			
	<b>MEB</b>	Bayonet	<input type="radio"/>	<input type="radio"/>
	<b>MEO</b>	Curved up	<input type="radio"/>	<input type="radio"/>
	<b>MEU</b>	Curved down	<input type="radio"/>	<input type="radio"/>

TIPS	TIP I.D.	STYLE	5mm X 30cm	5mm X 36cm
<b>DISSECTORS</b>				
	Cuschiери dissector			
	<b>KJB</b>	Bayonet	<input type="radio"/>	<input type="radio"/>
	<b>KJU</b>	Curved down	<input type="radio"/>	<input type="radio"/>
<b>SCISSORS</b>				
	Cuschiери serrated scissors			
	<b>MTB</b>	Bayonet	<input type="radio"/>	<input type="radio"/>
	<b>MTO</b>	Curved up	<input type="radio"/>	<input type="radio"/>
	<b>MTU</b>	Curved down	<input type="radio"/>	<input type="radio"/>
<b>BIOPSY &amp; PUNCH</b>				
	Cuschiери biopsy forceps			
	<b>DBB</b>	Bayonet	<input type="radio"/>	<input type="radio"/>
	<b>DBO</b>	Curved up	<input type="radio"/>	<input type="radio"/>
	<b>DBU</b>	Curved down	<input type="radio"/>	<input type="radio"/>

















# CLICK & go ARTICULATING

Articulating instruments make it possible for the surgeon to easily explore difficult to reach areas. A rotating, easy-motion wheel allows for continuous articulation up to 45°.

- 45° articulation
- 360° rotation
- Replaceable working inserts
- Easy dismantling for cleaning and sterilization
- Electrocautery capabilities
- Autoclavable

TIPS	KIT NUMBER	TIP INSERT NUMBER	OUTER SHEATH NUMBER	INNER SHEATH NUMBER	HANDLE NUMBER
<b>GRASPERS</b>					
	Grasping forceps, straight				
	<b>33521 AWL</b>	<b>33520 AWL</b>	<b>33543</b>	<b>33543 AL</b>	<b>33121</b>
	Duval grasping forceps				
	<b>33521 DUL</b>	<b>33520 DUL</b>	<b>33543</b>	<b>33543 AL</b>	<b>33121</b>
<b>SCISSORS</b>					
	Metzenbaum curved scissors, 12mm				
	<b>34521 MSL</b>	<b>34520 MSL</b>	<b>33543</b>	<b>33543 AL</b>	<b>33121</b>
<b>BIOPSY</b>					
	Biopsy forceps				
	<b>34521 BFL</b>	<b>34520 BFL</b>	<b>33543</b>	<b>33543 AL</b>	<b>33121</b>

















# CLICK & go ORDERING

				33121	33125	33126	33131	33132
								
				INSULATED KITS			NON-INSULATED KITS	
TIP I.D. (sorted alphabetically)	Dia.	Lgth	Insert†					
<b>A</b> 	3.5mm	30cm	30340 A	33221 A 33321 A 33421 A	33225 A 33325 A 33425 A	33226 A 33326 A 33426 A	30341 A	33232 A 33332 A 33432 A
	5mm	30cm	33210 A				33231 A	
	5mm	36cm	33310 A				33331 A	
	5mm	43cm	33410 A				33431 A	
<b>AF</b> 	5mm	30cm	33210 AF	33221 AF 33321 AF 33421 AF	33225 AF 33325 AF 33425 AF	33226 AF 33326 AF 33426 AF	33231 AF	33232 AF 33332 AF 33432 AF
	5mm	36cm	33310 AF				33331 AF	
	5mm	43cm	33410 AF				33431 AF	
<b>AG</b> 	5mm	30cm	33210 AG	33221 AG 33321 AG 33421 AG	33225 AG 33325 AG 33425 AG	33226 AG 33326 AG 33426 AG	33231 AG	33232 AG 33332 AG 33432 AG
	5mm	36cm	33310 AG				33331 AG	
	5mm	43cm	33410 AG				33431 AG	
<b>AR</b> 	5mm	30cm	33210 AR	33221 AR 33321 AR 33421 AR	33225 AR 33325 AR 33425 AR	33226 AR 33326 AR 33426 AR	33231 AR	33232 AR 33332 AR 33432 AR
	5mm	36cm	33310 AR				33331 AR	
	5mm	43cm	33410 AR				33431 AR	
<b>AS</b> 	2mm	20cm	30240 AS				30241 AS	
	3mm	20cm	30340 AS				30341 AS	
<b>AV</b> 	5mm	30cm	33210 AV				33231 AV	33232 AV 33332 AV 33432 AV
	5mm	36cm	33310 AV				33331 AV	
	5mm	43cm	33410 AV				33431 AV	
<b>BA</b>  Standard	5mm	30cm	33210 BA	33221 BA 33321 BA 33421 BA	33225 BA 33325 BA 33425 BA	33226 BA 33326 BA 33426 BA	33231 BA	33232 BA 33332 BA 33432 BA
	5mm	36cm	33310 BA				33331 BA	
	5mm	43cm	33410 BA				33431 BA	
<b>BAB</b>  Bayonet	5mm	30cm	35210 BAB		35261 BAB* 35361 BAB*			
	5mm	36cm	35310 BAB					
<b>BAO</b>  Curved up	5mm	30cm	35210 BAO		35261 BAO* 35361 BAO*			
	5mm	36cm	35310 BAO					
<b>BAU</b>  Curved Down	5mm	30cm	35210 BAU		35261 BAU* 35361 BAU*			
	5mm	36cm	35310 BAU					
<b>BB</b>  Standard	5mm	30cm	33210 BB	33221 BB 33321 BB 33421 BB	33225 BB 33325 BB 33425 BB	33226 BB 33326 BB 33426 BB	33231 BB	33232 BB 33332 BB 33432 BB
	5mm	36cm	33310 BB				33331 BB	
	5mm	43cm	33410 BB				33431 BB	

\* use Non-rotating 33125 P Metal, Insulated Handle.

† 2mm, 3mm and 3.5mm Clicklines include a single piece Insert/Outer Sheath assembly.



















				33121	33125	33126	33131	33132
					 Metal			
TIP I.D. (sorted alphabetically)	Dia.	Lgth	Insert†	INSULATED KITS			NON-INSULATED KITS	
<b>BBO</b>  Curved Up	5mm 5mm	30cm 36cm	35210 BBO 35310 BBO		35261 BBO* 35361 BBO*			
<b>BBU</b>  Curved Down	5mm 5mm	30cm 36cm	35210 BBU 35310 BBU		35261 BBU* 35361 BBU*			
<b>BC</b> 	10mm	36cm	33510 BC				33531 BC	33532 BC
<b>BL</b> 	10mm	36cm	33510 BL				33531 BL	33532 BL
<b>BLS</b> 	10mm	36cm	33510 BLS				33531 BLS	33532 BLS
<b>C</b> 	5mm 5mm 5mm	30cm 36cm 43cm	33210 C 33310 C 33410 C	33221 C 33321 C 33421 C	33225 C 33325 C 33425 C	33226 C 33326 C 33426 C	33231 C 33331 C 33431 C	33232 C 33332 C 33432 C
<b>CB</b> 	10mm	36cm	33510 CB				33531 CB	33532 CB
<b>CC</b> 	5mm 5mm 5mm	30cm 36cm 43cm	33210 CC 33310 CC 33410 CC	33221 CC 33321 CC 33421 CC	33225 CC 33325 CC 33425 CC	33226 CC 33326 CC 33426 CC	33231 CC 33331 CC 33431 CC	33232 CC 33332 CC 33432 CC
<b>CG</b> 	5mm 5mm 5mm	30cm 36cm 43cm	33210 CG 33310 CG 33410 CG	33221 CG 33321 CG 33421 CG	33225 CG 33325 CG 33425 CG	33226 CG 33326 CG 33426 CG	33231 CG 33331 CG 33431 CG	33232 CG 33332 CG 33432 CG
<b>D</b> 	5mm 5mm 5mm	30cm 36cm 43cm	33210 D 33310 D 33410 D	33221 D 33321 D 33421 D	33325 D 33325 D 33425 D	33226 D 33326 D 33426 D	33231 D 33331 D 33431 D	33232 D 33332 D 33432 D
<b>DB</b> 	3.5mm 5mm 5mm 5mm	30cm 30cm 36cm 43cm	30710 DB 34210 DB 34310 DB 34410 DB	30721 DB 34221 DB 34321 DB 34421 DB	34225 DB 34325 DB 34425 DB		34231 DB 34331 DB 34431 DB	

\* use Non-rotating 33125 P Metal, Insulated Handle.

† 2mm, 3mm and 3.5mm Clicklines include a single piece Insert/Outer Sheath assembly.

# CLICK & go ORDERING

















				33121	33125	33126	33131	33132
					 Metal			
TIP I.D. (sorted alphabetically)	Dia.	Lgth	Insert†	INSULATED KITS			NON-INSULATED KITS	
<b>DBB</b>  Bayonet	5mm 5mm	30cm 36cm						
<b>DBO</b>  Curved Up	5mm 5mm	30cm 36cm	35210 DBO 35310 DBO		35261 DBO* 35361 DBO*			
<b>DBS</b>  Curved Down	2mm 3mm	20cm 20cm	30210 DBS 30310 DBS	30221 DBS 30321 DBS				
<b>DBU</b>  Curved Down	5mm 5mm	30cm 36cm	35210 DBU 35310 DBU		35261 DBU* 35361 DBU*			
<b>DF</b> 	5mm 5mm 5mm	30cm 36cm 43cm	33210 DF 33310 DF 33410 DF	33221 DF 33321 DF 33421 DF	33225 DF 33325 DF 33425 DF	33226 DF 33326 DF 33426 DF	33231 DF 33331 DF 33431 DF	33232 DF 33332 DF 33432 DF
<b>DH</b> 	3.5mm	30cm	30710 DH	30721 DH				
<b>DHS</b> 	2mm 3mm	20cm 20cm	30210 DHS 30310 DHS	30221 DHS 30321 DHS				
<b>DK</b> 	5mm 5mm 5mm	30cm 36cm 43cm	33210 DK 33310 DK 33410 DK		33225DK 33325DK 33425DK	33226DK 33326DK 33426DK		
<b>DN</b> 	5mm 5mm 5mm	30cm 36cm 43cm	33210 DN 33310 DN 33410 DN	33221 DN 33321 DN 33421 DN	33225 DN 33325 DN 33425 DN	33226 DN 33326 DN 33426 DN	33231 DN 33331 DN 33431 DN	33232 DN 33332 DN 33432 DN
<b>DS</b> 	5mm 5mm 5mm	30cm 36cm 43cm	34210 DS 34310 DS 34410 DS	34221 DS 34321 DS 34421 DS	34225 DS 34325 DS 34425 DS		34231 DS 34331 DS 34431 DS	
<b>DU</b> 	10mm	36cm	33510 DU				33531 DU	33532 DU

\* use Non-rotating 33125 P Metal, Insulated Handle.

† 2mm, 3mm and 3.5mm Clicklines include a single piece Insert/Outer Sheath assembly.


















# CLICK & go ORDERING

TIP I.D. (sorted alphabetically)	Dia.	Lgth	Insert†	33121	33125	33126	33131	33132
								
				INSULATED KITS			NON-INSULATED KITS	
<b>DY</b> 	5mm 5mm 5mm	30cm 36cm 43cm	33210 DY 33310 DY 33410 DY	33221 DY 33321 DY 33421 DY	33225 DY 33325 DY 33425 DY	33226 DY 33326 DY 33426 DY	33231 DY 33331 DY 33431 DY	33232 DY 33332 DY 33432 DY
<b>DYD</b> 	5mm 5mm 5mm	30cm 36cm 43cm	33210 DYD 33310 DYD 33410 DYD	33221 DYD 33321 DYD 33421 DYD	33225 DYD 33325 DYD 33425 DYD	33226 DYD 33326 DYD 33426 DYD	33231 DYD 33331 DYD 33431 DYD	33232 DYD 33332 DYD 33432 DYD
<b>DZ</b> 	5mm 5mm 5mm	30cm 36cm 43cm	34210 DZ 34310 DZ 34410 DZ	34221 DZ 34321 DZ 34421 DZ	34225 DZ 34325 DZ 34425 DZ		34231 DZ 34331 DZ 34431 DZ	
<b>EH</b> 	3.5mm 5mm 5mm 5mm	30cm 30cm 36cm 43cm	30710 EH 34210 EH 34310 EH 34410 EH	30721 EH 34221 EH 34321 EH 34421 EH	34225 EH 34325 EH 34425 EH		34231 EH 34331 EH 34431 EH	
<b>EHS</b> 	3mm	20cm	30310 EHS	30321 EHS				
<b>EK</b> 	5mm 5mm 5mm	30cm 36cm 43cm	34210 EK 34310 EK 34410 EK	34221 EK 34321 EK 34421 EK	34225 EK 34325 EK 34425 EK		34231 EK 34331 EK 34431 EK	
<b>ES</b> 	5mm 5mm 5mm	30cm 36cm 43cm	34210 ES 34310 ES 34410 ES	34221 ES 34321 ES 34421 ES	34225 ES 34325 ES 34425 ES		34231 ES 34331 ES 34431 ES	
<b>FA</b> 	3mm	30cm	30340 FA				30331 FA**	
<b>FB</b> 	3mm	30cm	30340 FB				30331 FB**	
<b>FG</b> 	3.5mm 5mm 5mm 5mm	30cm 30cm 36cm 43cm	30710 FG 33210 FG 33310 FG 33410 FG	30721 FG 33221 FG 33321 FG 33421 FG	33225 FG 33325 FG 33425 FG	33226 FG 33326 FG 33426 FG	33231 FG 33331 FG 33431 FG	33232 FG 33332 FG 33432 FG
<b>FGS</b> 	3mm	20cm	30310 FGS	30321 FGS				

\*\* uses 30131 Mini Handle.

† 2mm, 3mm and 3.5mm Clicklines include a single piece Insert/Outer Sheath assembly.

















# CLICK & go ORDERING

				33121	33125	33126	33131	33132
								
				INSULATED KITS			NON-INSULATED KITS	
TIP I.D. (sorted alphabetically)	Dia.	Lgth	Insert†					
<b>FK</b> 	10mm	36cm	33510 FK				33531 FK	33532 FK
<b>FL</b> 	3mm	30cm	30340 FL				30331 FL**	
<b>FM</b> 	5mm 5mm 5mm	30cm 36cm 43cm	33210 FM 33310 FM 33410 FM	33221 FM 33321 FM 33421 FM	33225 FM 33325 FM 33425 FM	33226 FM 33326 FM 33426 FM	33231 FM 33331 FM 33431 FM	33232 FM 33332 FM 33432 FM
<b>FR</b> 	3mm	30cm	30340 FR				30331 FR**	
<b>FS</b> 	10mm	36cm	33510 FS				33531 FS	33532 FS
<b>G</b> 	5mm 5mm 5mm	30cm 36cm 43cm	33210 G 33310 G 33410 G	33221 G 33321 G 33421 G	33225 G 33325 G 33425 G	33226 G 33326 G 33426 G	33231 G 33331 G 33431 G	33232 G 33332 G 33432 G
<b>GF</b> 	10mm	36cm	33510 GF				33531 GF	33532 GF
<b>GS</b> 	10mm	36cm	34510 GS				34531 GS	
<b>HM</b> 	5mm 5mm 5mm	30cm 36cm 43cm	33210 HM 33310 HM 33410 HM	33221 HM 33321 HM 33421 HM	33225 HM 33325 HM 33425 HM	33226 HM 33326 HM 33426 HM	33231 HM 33331 HM 33431 HM	33232 HM 33332 HM 33432 HM
<b>J</b> 	5mm 5mm 5mm	30cm 36cm 43cm	33210 J 33310 J 33410 J	33221 J 33321 J 33421 J	33225 J 33325 J 33425 J	33226 J 33326 J 33426 J	33231 J 33331 J 33431 J	33232 J 33332 J 33432 J
<b>K</b> 	3.5mm 5mm 5mm 5mm	30cm 30cm 36cm 43cm	30710 K 33210 K 33310 K 33410 K	30721 K 33221 K 33321 K 33421 K	33225 K 33325 K 33425 K	33226 K 33326 K 33426 K	33231 K 33331 K 33431 K	33232 K 33332 K 33432 K

\*\* uses 30131 Mini Handle.

† 2mm, 3mm and 3.5mm Clicklines include a single piece Insert/Outer Sheath assembly.



















TIP I.D. (sorted alphabetically)	Dia.	Lgth	Insert†		 Metal			
				INSULATED KITS			NON-INSULATED KITS	
 KJ	5mm	30cm	33210 KJ	33221 KJ	33225 KJ	33226 KJ	33231 KJ	33232 KJ
	5mm	36cm	33310 KJ	33321 KJ	33325 KJ	33326 KJ	33331 KJ	33332 KJ
	5mm	43cm	33410 KJ	33421 KJ	33425 KJ	33426 KJ	33431 KJ	33432 KJ
 KJB  Bayonet	5mm	30cm	35210 KJB		35261 KJB*			
	5mm	36cm	35310 KJB		35361 KJB*			
 KJS	2mm	20cm	30210 KJS	30221 KJS				
 KJU  Curved Down	5mm	30cm	35210 KJU		35261 KJU*			
	5mm	36cm	35310 KJU		35361 KJU*			
 KS	3mm	20cm	30310 KS	30321 KS				
 KW	5mm	30cm	33210 KW	33221 KW	33225 KW	33226 KW	33231 KW	33232 KW
	5mm	36cm	33310 KW	33321 KW	33325 KW	33326 KW	33331 KW	33332 KW
	5mm	43cm	33410 KW	33421 KW	33425 KW	33426 KW	33431 KW	33432 KW
 LF	5mm	30cm	33210 LF	33221 LF	33225 LF	33226 LF	33231 LF	33232 LF
	5mm	36cm	33310 LF	33321 LF	33325 LF	33326 LF	33331 LF	33332 LF
	5mm	43cm	33410 LF	33421 LF	33425 LF	33426 LF	33431 LF	33432 LF
 LFD	5mm	30cm	33210 LFD	33221 LFD	33225 LFD	33226 LFD	33231 LFD	33232 LFD
	5mm	36cm	33310 LFD	33321 LFD	33325 LFD	33326 LFD	33331 LFD	33332 LFD
	5mm	43cm	33410 LFD	33421 LFD	33425 LFD	33426 LFD	33431 LFD	33432 LFD
 M	5mm	30cm	33210 M	33221 M	33225 M	33226 M	33231 M	33232 M
	5mm	36cm	33310 M	33321 M	33325 M	33326 M	33331 M	33332 M
	5mm	43cm	33410 M	33421 M	33425 M	33426 M	33431 M	33432 M
 MA	5mm	30cm	34210 MA	34221 MA	34225 MA		34231 MA	
	5mm	36cm	34310 MA	34321 MA	34325 MA		34331 MA	
	5mm	43cm	34410 MA	34421 MA	34425 MA		34431 MA	
 MA	5mm	30cm	33210 MA	33221 MA	33225 MA	33226 MA	33231 MA	33232 MA
	5mm	36cm	33310 MA	33321 MA	33325 MA	33326 MA	33331 MA	33332 MA
	5mm	43cm	33410 MA	33421 MA	33425 MA	33426 MA	33431 MA	33432 MA

\* use Non-rotating 33125 P Metal, Insulated Handle.

† 2mm, 3mm and 3.5mm Clicklines include a single piece Insert/Outer Sheath assembly.

















# CLICK & go ORDERING

				33121	33125	33126	33131	33132
								
				INSULATED KITS			NON-INSULATED KITS	
TIP I.D. (sorted alphabetically)	Dia.	Lgth	Insert†					
MAS	3mm	20cm	30310 MAS	30321 MAS				
								
MB	5mm	30cm	34210 MB	34221 MB	34225 MB		34231 MB	
	5mm	36cm	34310 MB	34321 MB	34325 MB		34331 MB	
	5mm	43cm	34410 MB	34421 MB	34425 MB		34431 MB	
MB	5mm	30cm	35210 MB		35261 MB*			
	5mm	36cm	35310 MB		35361 MB*			
Bayonet								
MD	5mm	30cm	34210 MD	34221 MD	34225 MD		34231 MD	
	5mm	36cm	34310 MD	34321 MD	34325 MD		34331 MD	
	5mm	43cm	34410 MD	34421 MD	34425 MD		34431 MD	
MD	3.5mm	30cm	30710 MD	30721 MD				
	5mm	30cm	33210 MD	33221 MD	33225 MD	33226 MD	33231 MD	33232 MD
	5mm	36cm	33310 MD	33321 MD	33325 MD	33326 MD	33331 MD	33332 MD
	5mm	43cm	33410 MD	33421 MD	33425 MD	33426 MD	33431 MD	33432 MD
MDS	3mm	20cm	30310 MDS	30321 MDS				
								
ME	5mm	30cm	33210 ME	33221 ME	33225 ME	33226 ME	33231 ME	33232 ME
	5mm	36cm	33310 ME	33321 ME	33325 ME	33326 ME	33331 ME	33332 ME
Standard	5mm	43cm	33410 ME	33421 ME	33425 ME	33426 ME	33431 ME	33432 ME
MEB	5mm	30cm	35210 MEB		35261 MEB*			
	5mm	36cm	35310 MEB		35361 MEB*			
Bayonet								
MEM	5mm	30cm	33210 MEM	33221 MEM	33225 MEM	33226 MEM	33231 MEM	33232 MEM
	5mm	36cm	33310 MEM	33321 MEM	33325 MEM	33326 MEM	33331 MEM	33332 MEM
	5mm	43cm	33410 MEM	33421 MEM	33425 MEM	33426 MEM	33431 MEM	33432 MEM
MEO	5mm	30cm	35210 MEO		35261 MEO*			
	5mm	36cm	35310 MEO		35361 MEO*			
Curved Up								
MES	3mm	20cm	30310 MES	30321 MES				
								

\* use Non-rotating 33125 P Metal, Insulated Handle.

† 2mm, 3mm and 3.5mm Clicklines include a single piece Insert/Outer Sheath assembly.



















TIP I.D. (sorted alphabetically)	Dia.	Lgth	Insert†	33121	33125	33126	33131	33132
					 Metal			
<b>MEU</b>  Curved Down	5mm 5mm	30cm 36cm	35210 MEU 35310 MEU	INSULATED KITS			NON-INSULATED KITS	
					35261 MEU* 35361 MEU*			
<b>MF</b> 	5mm 5mm 5mm	30cm 36cm 43cm	33210 MF 33310 MF 33410 MF	33221 MF 33321 MF 33421 MF	33225 MF 33325 MF 33425 MF	33226 MF 33326 MF 33426 MF	33231 MF 33331 MF 33431 MF	33232 MF 33332 MF 33432 MF
<b>MG</b> 	5mm 5mm 5mm	30cm 36cm 43cm	33210 MG 33310 MG 33410 MG	33221 MG 33321 MG 33421 MG	33225 MG 33325 MG 33425 MG	33226 MG 33326 MG 33426 MG	33231 MG 33331 MG 33431 MG	33232 MG 33332 MG 33432 MG
<b>MGS</b> 	3mm	20cm	30310 MGS	30321 MGS				
<b>MH</b> 	5mm 5mm 5mm	30cm 36cm 43cm	33210 MH 33310 MH 33410 MH	33221 MH 33321 MH 33421 MH	33225 MH 33325 MH 33425 MH	33226 MH 33326 MH 33426 MH	33231 MH 33331 MH 33431 MH	33232 MH 33332 MH 33432 MH
<b>MHS</b> 	3mm	20cm	30310 MHS	30321 MHS				
<b>MKM</b> 	5mm 5mm 5mm	30cm 36cm 43cm	33210 MKM 33310 MKM 33410 MKM	33221 MKM 33321 MKM 33421 MKM	33225 MKM 33325 MKM 33425 MKM	33226 MKM 33326 MKM 33426 MKM	33231 MKM 33331 MKM 33431 MKM	33232 MKM 33332 MKM 33432 MKM
<b>ML</b> 	5mm 5mm 5mm	30cm 36cm 43cm	33210 ML 33310 ML 33410 ML	33221 ML 33321 ML 33421 ML	33225 ML 33325 ML 33425 ML	33226 ML 33326 ML 33426 ML	33231 ML 33331 ML 33431 ML	33232 ML 33332 ML 33432 ML
<b>MLL</b> 	10mm	36cm	33510 MLL				33531 MLL	33532 MLL
<b>MLM</b> 	5mm 5mm 5mm	30cm 36cm 43cm	33210 MLM 33310 MLM 33410 MLM	33221 MLM 33321 MLM 33421 MLM	33225 MLM 33325 MLM 33425 MLM	33226 MLM 33326 MLM 33426 MLM	33231 MLM 33331 MLM 33431 MLM	33232 MLM 33332 MLM 33432 MLM
<b>MLS</b> 	3mm	20cm	30310 MLS	30321 MLS				

\* use Non-rotating 33125 P Metal, Insulated Handle.

† 2mm, 3mm and 3.5mm Clicklines include a single piece Insert/Outer Sheath assembly.

















# CLICK & go ORDERING

				33121	33125	33126	33131	33132
								
TIP I.D. (sorted alphabetically)	Dia.	Lgth	Insert†	INSULATED KITS			NON-INSULATED KITS	
 <b>MNM</b>	5mm	30cm	33210 MNM	33221 MNM	33225 MNM	33226 MNM	33231 MNM	33232 MNM
	5mm	36cm	33310 MNM	33321 MNM	33325 MNM	33326 MNM	33331 MNM	33332 MNM
	5mm	43cm	33410 MNM	33421 MNM	33425 MNM	33426 MNM	33431 MNM	33432 MNM
 <b>MO</b> Curved Up	5mm	30cm	35210 MO		35261 MO*			
	5mm	36cm	35310 MO		35361 MO*			
 <b>MP</b>	5mm	30cm	33210 MP	33221 MP	33225 MP	33226 MP	33231 MP	33232 MP
	5mm	36cm	33310 MP	33321 MP	33325 MP	33326 MP	33331 MP	33332 MP
	5mm	43cm	33410 MP	33421 MP	33425 MP	33426 MP	33431 MP	33432 MP
 <b>MR</b>	5mm	30cm	33210 MR	33221 MR	33225 MR	33226 MR	33231 MR	33232 MR
	5mm	36cm	33310 MR	33321 MR	33325 MR	33326 MR	33331 MR	33332 MR
	5mm	43cm	33410 MR	33421 MR	33425 MR	33426 MR	33431 MR	33432 MR
 <b>MS</b>	5mm	30cm	33210 MS	33221 MS	33225 MS	33226 MS	33231 MS	33232 MS
	5mm	36cm	33310 MS	33321 MS	33325 MS	33326 MS	33331 MS	33332 MS
	5mm	43cm	33410 MS	33421 MS	33425 MS	33426 MS	33431 MS	33432 MS
 <b>MS</b>	5mm	30cm	34210 MS	34221 MS	34225 MS		34231 MS	
	5mm	36cm	34310 MS	34321 MS	34325 MS		34331 MS	
	5mm	43cm	34410 MS	34421 MS	34425 MS		34431 MS	
 <b>MSP</b>	2mm	20cm	30210 MSP	30221 MSP				
	3mm	20cm	30310 MSP	30321 MSP				
 <b>MSS</b>		20cm	30210 MSS	30221 MSS				
 <b>MT</b>	5mm	30cm	33210 MT	33221 MT	33225 MT	33226 MT	33231 MT	33232 MT
	5mm	36cm	33310 MT	33321 MT	33325 MT	33326 MT	33331 MT	33332 MT
	5mm	43cm	33410 MT	33421 MT	33425 MT	33426 MT	33431 MT	33432 MT
 <b>MT</b>	5mm	30cm	34210 MT	34221 MT	34225 MT		34231 MT	
	5mm	36cm	34310 MT	34321 MT	34325 MT		34331 MT	
	5mm	43cm	34410 MT	34421 MT	34425 MT		34431 MT	
 <b>MTB</b> Bayonet	5mm	30cm	35210 MTB		35261 MTB*			
	5mm	36cm	35310 MTB		35361 MTB*			

\* use Non-rotating 33125 P Metal, Insulated Handle.

† 2mm, 3mm and 3.5mm Clicklines include a single piece Insert/Outer Sheath assembly.



















				33121	33125	33126	33131	33132
					 Metal			
TIP I.D. (sorted alphabetically)	Dia.	Lgth	Insert†	INSULATED KITS			NON-INSULATED KITS	
<b>MTO</b>  Curved Up	5mm 5mm	30cm 36cm	35210 MTO 35310 MTO		35261 MTO* 35361 MTO*			
<b>MTU</b>  Curved Down	5mm 5mm	30cm 36cm	35210 MTU 35310 MTU		35261 MTU* 35361 MTU*			
<b>MU</b>  Curved Down	5mm 5mm 5mm	30cm 36cm 43cm	33210 MU 33310 MU 33410 MU	33221 MU 33321 MU 33421 MU	33225 MU 33325 MU 33425 MU	33226 MU 33326 MU 33426 MU	33231 MU 33331 MU 33431 MU	33232 MU 33332 MU 33432 MU
<b>MU</b>  Curved Down	5mm 5mm	30cm 36cm	35210 MU 35310 MU		35261 MU* 35361 MU*			
<b>MV</b>  Curved Down	5mm 5mm 5mm	30cm 36cm 43cm	33210 MV 33310 MV 33410 MV	33221 MV 33321 MV 33421 MV	33225 MV 33325 MV 33425 MV	33226 MV 33326 MV 33426 MV	33231 MV 33331 MV 33431 MV	33232 MV 33332 MV 33432 MV
<b>MW</b>  Curved Down	5mm 5mm 5mm	30cm 36cm 43cm	33210 MW 33310 MW 33410 MW	33221 MW 33321 MW 33421 MW	33225 MW 33325 MW 33425 MW	33226 MW 33326 MW 33426 MW	33231 MW 33331 MW 33431 MW	33232 MW 33332 MW 33432 MW
<b>MW</b>  Curved Down	3.5mm 5mm 5mm 5mm	30cm 30cm 36cm 43cm	30710 MW 34210 MW 34310 MW 34410 MW	30721 MW 34221 MW 34321 MW 34421 MW	34225 MW 34325 MW 34425 MW		34231 MW 34331 MW 34431 MW	
<b>MWS</b>  Curved Down	3mm	20cm	30310 MWS	30321 MWS				
<b>NS</b>  Curved Down	3mm	20cm	30310 NS	30321 NS				
<b>ON</b>  Curved Down	3.5mm 5mm 5mm 5mm	30cm 30cm 36cm 43cm	30710 ON 33210 ON 33310 ON 33410 ON	30721 ON 33221 ON 33321 ON 33421 ON	33225 ON 33325 ON 33425 ON	33226 ON 33326 ON 33426 ON	33231 ON 33331 ON 33431 ON	33232 ON 33332 ON 33432 ON
<b>ONS</b>  Curved Down	3mm	20cm	30310 ONS	30321 ONS				

\* use Non-rotating 33125 P Metal, Insulated Handle.

† 2mm, 3mm and 3.5mm Clicklines include a single piece Insert/Outer Sheath assembly.

















# CLICK & go ORDERING

				33121	33125	33126	33131	33132
					 Metal			
TIP I.D. (sorted alphabetically)	Dia.	Lgth	Insert†	INSULATED KITS			NON-INSULATED KITS	
PG 	10mm	36cm	33510 PG				33531 PG	
PR 	10mm	36cm	33510 PR				33531 PR	33532PR
R 	5mm 5mm 5mm	30cm 36cm 43cm	33210 R 33310 R 33410 R	33221 R 33321 R 33421 R	33225 R 33325 R 33425 R	33226 R 33326 R 33426 R	33231 R 33331 R 33431 R	33232 R 33332 R 33432 R
RA 	5mm 5mm 5mm	30cm 36cm 43cm	33210 RA 33310 RA 33410 RA	33221 RA 33321 RA 33421 RA	33225 RA 33325 RA 33425 RA	33226 RA 33326 RA 33426 RA	33231 RA 33331 RA 33431 RA	33232 RA 33332 RA 33432 RA
RG 	10mm	36cm	33510 RG				33531 RG	33532 RG
RM 	5mm 5mm 5mm	30cm 36cm 43cm	33210 RM 33310 RM 33410 RM	33221 RM 33321 RM 33421 RM	33225 RM 33325 RM 33425 RM	33226 RM 33326 RM 33426 RM		
RS 	3mm	20cm	30310 RS	30321 RS				
S 	5mm 5mm	30cm 36cm	34210 S 34310 S	34221 S 34321 S	34225 S 34325 S		34231 S 34331 S	
SA 	3mm	30cm	30340 SA				30331 SA**	
SB 	3mm	30cm	30340 SB				30331 SB**	
SBS 	3mm	20cm	30340 SBS				30331 SBS	

\*\* uses 30131 Mini Handle.

† 2mm, 3mm and 3.5mm Clicklines include a single piece Insert/Outer Sheath assembly.




















				33121	33125	33126	33131	33132
								
				Metal				
TIP I.D.	Dia.	Lgth	Insert†	INSULATED KITS			NON-INSULATED KITS	
(sorted alphabetically)								
	5mm	30cm	33210 SD	33221 SD	33225 SD	33226 SD	33231 SD	33232 SD
	5mm	36cm	33310 SD	33321 SD	33325 SD	33326 SD	33331 SD	33332 SD
	5mm	43cm	33410 SD	33421 SD	33425 SD	33426 SD	33431 SD	33432 SD
	5mm	30cm	34210 SF	34221 SF	34225 SF		34231 SF	
	5mm	36cm	34310 SF	34321 SF	34325 SF		34331 SF	
	10mm	36cm	33510 SG				33531 SG	33532 SG
	5mm	30cm	34210 SH	34221 SH	34225 SH	34226 SH		
	5mm	36cm	34310 SH	34321 SH	34325 SH	34326 SH		
	5mm	43cm	34410 SH	34421 SH	34425 SH	34426 SH		
	5mm	30cm	34210 SK	34221 SK	34225 SK		34231 SK	
	5mm	36cm	34310 SK	34321 SK	34325 SK		34331 SK	
	5mm	30cm	34210 SL	34221 SL	34225 SL		34231 SL	
	5mm	36cm	34310 SL	34321 SL	34325 SL		34331 SL	
	5mm	30cm	33210 SN	33221 SN	33225 SN	33226 SN	33231 SN	33232 SN
	5mm	36cm	33310 SN	33321 SN	33325 SN	33326 SN	33331 SN	33332 SN
	5mm	43cm	33410 SN	33421 SN	33425 SN	33426 SN	33431 SN	33432 SN
	5mm	30cm	33210 SW	33221 SW	33225 SW	33226 SW	33231 SW	33232 SW
	5mm	36cm	33310 SW	33321 SW	33325 SW	33326 SW	33331 SW	33332 SW
	5mm	43cm	33410 SW	33421 SW	33425 SW	33426 SW	33431 SW	33432 SW
	5mm	30cm	33210 TF	33221 TF	33225 TF	33226 TF	33231 TF	33232 TF
	5mm	36cm	33310 TF	33321 TF	33325 TF	33326 TF	33331 TF	33332 TF
	5mm	43cm	33410 TF	33421 TF	33425 TF	33426 TF	33431 TF	33432 TF
	10mm	36cm	33510 TH				33531 TH	33532 TH
	3.5mm	30cm	30710 UL	30721 UL				
	5mm	30cm	33210 UL	33221 UL	33225 UL	33226 UL	33231 UL	33232 UL
	5mm	36cm	33310 UL	33321 UL	33325 UL	33326 UL	33331 UL	33332 UL
	5mm	43cm	33410 UL	33421 UL	33425 UL	33426 UL	33431 UL	33432 UL

\* use Non-rotating 33125 P Metal, Insulated Handle.

† 2mm, 3mm and 3.5mm Clicklines include a single piece Insert/Outer Sheath assembly.

# CLICK & go ORDERING

				33121	33125	33126	33131	33132
								
TIP I.D.	Dia.	Lgth	Insert†	INSULATED KITS			NON-INSULATED KITS	
(sorted alphabetically)								
<b>ULS</b>	3mm	20cm	30310 ULS	30321 ULS				
								
<b>UM</b>	5mm	30cm	33210 UM				33231 UM	33232 UM
	5mm	36cm	33310 UM				33331 UM	33332 UM
	5mm	43cm	33410 UM				33431 UM	33432 UM
	10mm	36cm	33510 UM				33531 UM	33532 UM
<b>UN</b>	10mm	36cm	33510 UN				33531 UN	33532 UN
								
<b>UO</b>	10mm	36cm	33510 UO				33531 UO	33532 UO
								
<b>VA</b>	5mm	30cm	33210 VA	33221 VA	33225 VA	33226 VA	33231 VA	33232 VA
	5mm	36cm	33310 VA	33321 VA	33325 VA	33326 VA	33331 VA	33332 VA
	5mm	43cm	33410 VA	33421 VA	33425 VA	33426 VA	33431 VA	33432 VA
<b>VDS</b>	3mm	20cm	30340 VDS				30341 VDS	
								
<b>VK</b>	5mm	30cm	33210 VK	33221 VK	33225 VK	33226 VK	33231 VK	33232 VK
	5mm	36cm	33310 VK	33321 VK	33325 VK	33326 VK	33331 VK	33332 VK
	5mm	43cm	33410 VK	33421 VK	33425 VK	33426 VK	33431 VK	33432 VK
<b>VT</b>	5mm	30cm	33210 VT	33221 VT	33225 VT	33226 VT	33231 VT	33232 VT
	5mm	36cm	33310 VT	33321 VT	33325 VT	33326 VT	33331 VT	33332 VT
	5mm	43cm	33410 VT	33421 VT	33425 VT	33426 VT	33431 VT	33432 VT
<b>VRS</b>	3mm	20cm	30340 VRS				30341 VRS	
								
<b>WS</b>	10mm	36cm	34510 WS				34531 WS	
								
<b>WT</b>	5mm	30cm	33210 WT	33221 WT	33225 WT	33226 WT	33231 WT	33232 WT
	5mm	36cm	33310 WT	33321 WT	33325 WT	33326 WT	33331 WT	33332 WT
	5mm	43cm	33410 WT	33421 WT	33425 WT	33426 WT	33431 WT	33432 WT
<b>WTD</b>	5mm	30cm	33210 WTD	33221 WTD	33225 WTD	33226 WTD	33231 WTD	33232 WTD
	5mm	36cm	33310 WTD	33321 WTD	33325 WTD	33326 WTD	33331 WTD	33332 WTD
	5mm	43cm	33410 WTD	33421 WTD	33425 WTD	33426 WTD	33431 WTD	33432 WTD

† 2mm, 3mm and 3.5mm Clicklines include a single piece Insert/Outer Sheath assembly.





# CLICK & go TELESCOPES

## STANDARD LENGTH 5MM TELESCOPES

<b>26046 AA</b>	Hopkins II telescope, 0°, enlarged view, 5mm x 29cm, autoclavable, improved image size, optical clarity and brightness
<b>26046 BA</b>	Hopkins II telescope, 30°, enlarged view, 5mm x 29cm, autoclavable, improved image size, optical clarity and brightness
<b>26046 FA</b>	Hopkins II telescope, 45°, enlarged view, 5mm x 29cm, autoclavable, improved image size, optical clarity and brightness
<b>26006 AA</b>	Hopkins II telescope, 0°, enlarged view, 5mm x 29cm, autoclavable
<b>26006 BA</b>	Hopkins II telescope, 30°, enlarged view, 5mm x 29cm, autoclavable
<b>26006 FA</b>	Hopkins II telescope, 45°, enlarged view, 5mm x 29cm, autoclavable

## STANDARD LENGTH 10MM TELESCOPES

<b>26003 AA</b>	Hopkins II telescope, 0°, enlarged view, 10mm x 31cm, autoclavable
<b>26003 BA</b>	Hopkins II telescope, 30°, enlarged view, 10mm x 31cm, autoclavable
<b>26003 FA</b>	Hopkins II telescope, 45°, enlarged view, 10mm x 31cm, autoclavable
<b>26003 AE</b>	EndoCAMEleon telescope, 10mm x 32cm, variable DOV from 0°-120°

## EXTENDED LENGTH 10MM TELESCOPES

<b>26003 AEA</b>	Hopkins II telescope, 0°, enlarged view, 10mm x 42cm, autoclavable
<b>26003 BEA</b>	Hopkins II telescope, 30°, enlarged view, 10mm x 42cm, autoclavable
<b>26003 FEA</b>	Hopkins II telescope, 45°, enlarged view, 10mm x 42cm, autoclavable
<b>26003 AEE</b>	EndoCAMEleon telescope, 10mm x 42cm, variable DOV from 0°-120°

## MINI & PEDIATRIC LENGTH TELESCOPES

<b>26011 AA</b>	Hopkins II telescope, 0°, enlarged view, 5mm x 24cm, autoclavable
<b>26011 BA</b>	Hopkins II telescope, 30°, enlarged view, 5mm x 24cm, autoclavable
<b>26009 AA</b>	Hopkins II telescope, 0°, 4mm x 18cm, autoclavable
<b>26009 BA</b>	Hopkins II telescope, 30°, 4mm x 18cm, autoclavable
<b>26007 AA</b>	Hopkins II telescope, 0°, enlarged view, 3.3mm x 25cm, autoclavable
<b>26007 BA</b>	Hopkins II telescope, 30°, enlarged view, 3.3mm x 25cm, autoclavable
<b>26007 FA</b>	Hopkins II telescope, 45°, enlarged view, 3.3mm x 25cm, autoclavable
<b>26008 AA</b>	Hopkins II telescope, 0°, 2mm x 26cm, autoclavable
<b>26008 BA</b>	Hopkins II telescope, 30°, 2mm x 26cm, autoclavable
<b>11602 AA</b>	Miniature telescope, semi-rigid, 0°, 2mm x 20cm, autoclavable

# CLICK & go OPERATIVE TELESCOPES

<b>26034 AAK</b>	Hopkins telescope, 0°, with parallel eyepiece, 10mm x 23cm, autoclavable. Fiber optic light transmission incorporated, with 6mm instrument channel.	<b>26038 AA</b>	Hopkins telescope, 0°, with angled eyepiece, 10mm x 27cm, autoclavable. Fiber optic light transmission incorporated, with 6mm instrument channel.
<b>26036 AA</b>	Hopkins telescope, 0°, with parallel eyepiece, enlarged view, 10mm x 25.5cm, autoclavable. Fiber optic light transmission incorporated, with 5mm instrument channel.	<b>26075 AA</b>	Hopkins telescope, 6°, with parallel eyepiece, wide angle, 11mm x 27.5cm, autoclavable. Fiber optic light transmission incorporated, with 7.5mm instrument channel.
		<b>26075 AAS</b>	Hopkins telescope, 6°, with parallel eyepiece, wide angle, 11mm x 27.5cm, with standard insufflation stopcock, autoclavable. Fiber optic light transmission incorporated, with 7.5mm instrument channel.

# THE FUTURE OF MEDICINE STARTS WITH ONE.

KARL STORZ brings you the future of medicine with products that exceed your expectations. As your single source solution, we offer one of the broadest selections of rigid and flexible endoscopes, hand instruments, insufflators, suction and irrigation devices, video cameras, illumination sources and documentation solutions.

## OR1®

KARL STORZ brings you the integrated operating room where virtually every operating room component functions as a single, optimized system.

## IMAGE 1®

The Image 1® platform now includes Full HD endoscopic cameras for the best viewing experience. Our new 16:9 CCD chips capture HD images that are displayed on our new WideView™ 16:9 monitors with 1080p resolution—the highest standard in HDTV.

## PROTECTION 1®

And KARL STORZ works to add value at every stage of the customer relationship with comprehensive service and support, attractive financing options and asset protection through our Protection 1® Repair Program.



For more information on KARL STORZ products please contact us at 1-800-421-0837.

**STORZ**  
KARL STORZ—ENDOSCOPE  
THE DIAMOND STANDARD

KARL STORZ  
Endoscopy-America, Inc.

2151 E. Grand Ave  
El Segundo, CA 90245-5017

Phone 424 218 8100  
Toll Free 800 421 0837

Toll Free In Canada 800 268 4880  
[www.karlstorz.com](http://www.karlstorz.com)

[www.karlstorz.com](http://www.karlstorz.com)