



MiniLap® Percutaneous Surgical System


Less trauma. In its simplest form.

Introducing the MiniLap Percutaneous Surgical System from Teleflex. **Designed to help surgeons make minimally invasive surgery even less invasive**, the MiniLap System allows surgeons to insert the 2.4 mm shaft diameter instrument percutaneously through the skin using an integrated needle tip. The goal? Smaller incisions and less tissue trauma while maintaining the practices and techniques surgeons already use.

Less Trauma

- The 2.4 mm shaft allows for insertion directly through the skin, leaving a nearly unnoticeable scar
- Without the use of a trocar, the MiniLap System can provide the ability to add devices for triangulation and exposure
- The integrated needle tip allows percutaneous insertion without trocars
- In many cases, 2 or more trocars can be eliminated

Defect Comparison



Benefits




The Clinician
Uses the technique you already know




The Patient
Less trauma.

MiniLap® System With MiniGrip® Handle



MiniLap® System With ThumbGrip Handle



MiniLap System with MiniGrip Handle

ITEM NUMBER	PRODUCT DESCRIPTION	UNIT PER BOX
PGAC300	Alligator Grasper	6
PGCC300	Clutch Grasper	6
PGBK300	Babcock Grasper	6
PGBC300	Tong Grasper	6

MiniLap System with ThumbGrip Handle

ITEM NUMBER	PRODUCT DESCRIPTION	UNIT PER BOX
GBC250	Alligator Grasper	6
CLC250	Clutch Grasper	6
BCK250	Babcock Grasper	6
BLC250	Bowel Grasper	6
ECMC300	Curved Spatula Probe	6
ECMS300	Straight Spatula Probe	6
ECMH300	Hook Probe	6
ECMP300	Conical Probe	6

teleflex.com/minilap



Teleflex, the Teleflex logo, MiniGrip and MiniLap are trademarks or registered trademarks of Teleflex Incorporated or its affiliates, in the U.S. and/or other countries. Information in this document is not a substitute for the product Instructions for Use. The products in this document may not be available in all countries. Please contact your local representative. All data current at time of printing (07/2016). Subject to technical changes without further notice. © 2016 Teleflex Incorporated. All rights reserved. MC-000829 Rev. 2

Teleflex®

Teleflex · 3015 Carrington Mill Boulevard, Morrisville, NC 27560
Toll Free: 866 246 6990 | Phone: +1 919 544 8000
teleflex.com